


Cell cultivation and immunofluorescence staining with μ -Slide VI 0.4


In this protocol we describe a single example of cultivating HT-1080 cancer cells inside the μ -Slide VI^{0.4}. Subsequently, we stained the F-actin cytoskeleton with Alexa Fluor® 488 phalloidin and counterstained the nucleus with DAPI.

The protocol consists of four main steps:



1) Cultivation

- Unpack a μ -Slide VI^{0.4}, ibiTreat (80606) under sterile conditions and put it on a μ -Slide rack (80003). Apply 30 μ l of a 3×10^5 cells/ml HT-1080 cell suspension into each channel. Pipet directly into the channel as illustrated below or shown on our website.



- Cover reservoirs with the supplied lid.
- Put the slides with the rack into the incubator (37 °C; 5 % CO₂) and let cells attach (60 min). Afterwards fill both reservoirs with 60 μ l of cell-free medium.
- Incubate over night.

Application Note 09

2) Fixation of cells

- Aspirate medium from all reservoirs using a cell culture aspiration device. Wash cells with Dulbecco's PBS by slowly applying 200 μ l into one empty reservoir of each channel and aspirating from the opposite reservoir for each channel. Don't aspirate the entire channel volume.
- Fix cells with ~100 μ l of 3.7 % para-formaldehyde in PBS. After 10 min flush the liquid inside the channel by filling one well with 200 μ l PBS and removing the content of the reservoir from the other well; ensuring the channel is never dry.

Permeabilization and blocking

- Wash cells again with 200 μ l PBS as described above.
- Apply ~100 μ l of 0.1% Triton® X-100 (Fluka) in PBS for 3-5 min.
- Wash cells with PBS.
- Apply ~100 μ l of 1% BSA in PBS solution for 20 min.
- Wash cells with PBS.

3) Staining

- **Remove all liquid from the channel using your aspiration device. Don't let the channel dry.**
- Right after that, apply **25 μ l** of Alexa Fluor® 488 phalloidin (1 Unit + 500 μ l PBS + 1% BSA, Invitrogen Corp.) Incubate at room temperature for 20 min.
- Wash cells with PBS.
- Apply **25 μ l** DAPI (0.1 μ g/ml, Sigma-Aldrich) for 3-5 min.
- Wash cells with PBS and apply ibidi Mounting Medium until the channel is filled (ca. 50 μ l). ibidi Mounting Medium is glycerol based and contains DABCO for anti-fading*. The slide can be stored for approx. 4 weeks.

4) Imaging

- Observe cells under a fluorescence microscope with appropriate filter sets and optionally with immersion oil.



* As published: Lee J.H., Koh H., Kim M., Kim Y., Lee S.Y., Karess R.E., Lee S.-H., Shong M., Kim J.-M., Kim J. & Chung J.; Energy-dependent regulation of cell structure by AMP-activated protein kinase. Nature Aug **2007**; doi:10.1038/nature05828